

General 5-Hydroxyindoleacetic Acid (5-HIAA) ELISA Kit

Catalog No: #EK20544

Package Size: #EK20544-1 48T #EK20544-2 96T



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | General 5-Hydroxyindoleacetic Acid (5-HIAA) ELISA Kit |
| Applications | ELISA |
| Species Reactivity | General |

Application Details

Detect Range:6.25-400ng/mL

Sensitivity:2.93ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Detection Method:competitive inhibition

Note: This product is for in vitro research use only and is not intended for use in humans or animals.