

# Human Luteinizing hormone (LH) ELISA Kit

Catalog No: #EK20572



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human Luteinizing hormone (LH) ELISA Kit
Applications	ELISA
Species Reactivity	Human

## Application Details

Detect Range: 1 mIU/mL - 50 mIU/mL.

Sensitivity: 0.5 mIU/mL

Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50  $\mu$ L

Assay Time: 1-4.5h

Detection wavelength: 450 nm

Detection Method: sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.