

Recombinant Human Neuregulin 1-alpha EGF Domain

Catalog No: #AP60115




Package Size: #AP60115-1 10ug #AP60115-2 100ug #AP60115-3 500ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Neuregulin 1-alpha EGF Domain
Host Species	Escherichia coli.
Purification	> 97 % by SDS-PAGE and HPLC analyses.
Calculated MW	Approximately 7.4 kDa, a single non-glycosylated polypeptide chain containing 65 amino acids.
Target Sequence	SHLVKCAEKE KTFCVNGGEC FMVKDLSNPS RYLCKCQPGF TGARCTENVP MKVQNQEKAE ELYQK
Formulation	Lyophilized from a 0.2 µm filtered solution in 20 mM PB, pH 6.0, 150 mM NaCl.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. -□ A minimum of 12 months from date of receipt, when stored at ≤-20 °C as supplied. -□ 1 month, 2 to 8 °C under sterile conditions after reconstitution. -□ 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Images



Note: This product is for in vitro research use only and is not intended for use in humans or animals.