

Recombinant Murine Interleukin-21

Catalog No: #AP60163

Package Size: #AP60163-1 10ug #AP60163-2 100ug #AP60163-3 500ug


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Murine Interleukin-21
Host Species	Escherichia coli.
Purification	> 98 % by SDS-PAGE and HPLC analyses.
Calculated MW	Approximately 15.0 kDa, a single non-glycosylated polypeptide chain containing 129 amino acids.
Target Sequence	HKSSPQGPDR LLIRLRHLID IVEQLKIYEN DLDPELLSAP QDVKGHCHEA AFACFQKAKL KPSNPGNNKT FIIDLVAQLR RRLPARRGGK KQKHIACPS CDSYEK RTPK EFLERLKWLL QKMIHQHLS
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in 2 x PBS, pH 7.4.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. <ul style="list-style-type: none"> -□ A minimum of 12 months from date of receipt, when stored at &le;-20 °C as supplied. -□ 1 month, 2 to 8 °C under sterile conditions after reconstitution. -□ 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Images



Note: This product is for in vitro research use only and is not intended for use in humans or animals.