

Recombinant Murine LIGHT/TNFSF14

Catalog No: #AP60181


Package Size: #AP60181-1 5ug #AP60181-2 100ug #AP60181-3 500ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Murine LIGHT/TNFSF14
Host Species	Escherichia coli.
Purification	> 96 % by SDS-PAGE and HPLC analyses.
Calculated MW	Approximately 18.4 kDa, a single non-glycosylated polypeptide chain containing 168 amino acids.
Target Sequence	DGGKGSWEKL IQDQRSHQAN PAAHLTGANA SLIGIGGPLL WETRLGLAFL RGLTYHDGAL VTMEPGYYYYV YSKVQLSGVG CPQGLANGLP ITHGLYKRTS RYPKELELLV SRRSPCGRAN SSRVWWDSSF LGGVHLEAG EEVVVRVPGN RLVPRPDGTR SYFGAFMV
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH7.4.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. <input type="checkbox"/> A minimum of 12 months from date of receipt, when stored at ≤-20 °C as supplied. <input type="checkbox"/> 1 month, 2 to 8 °C under sterile conditions after reconstitution. <input type="checkbox"/> 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Images



Note: This product is for in vitro research use only and is not intended for use in humans or animals.