

Recombinant Murine Interferon-lambda2/Interleukin-28A

Catalog No: #AP60200


Package Size: #AP60200-1 5ug #AP60200-2 100ug #AP60200-3 500ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Murine Interferon-lambda2/Interleukin-28A
Host Species	E.coli
Purification	> 95 % by SDS-PAGE and HPLC analyses.
Calculated MW	Approximately 19.7 kDa, a single non-glycosylated polypeptide chain containing 174 amino acids.
Target Sequence	DPVPRATRLP VEAKDCHIAQ FKSLSPKELQ AFKKAKDAIE KRLLEKDMRC SSSLISRAWD LKQLVQERP KALQAEVALT LKVWENMTDS ALATILGQPL HTLSHIHSQL QTCTQLQATA EPKPPSRRLS RWLHRLQEAQ SKETPGCLED SVTSNLFRL TRDLKCVASG DQCV
Formulation	LyophilizedB fromB aB 0.2B umB filteredB concentratedB solutionB inB PBS.B
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. -□A minimum of 12 months from date of receipt, when stored at ≤-20 °C as supplied. -□1 month, 2 to 8 °C under sterile conditions after reconstitution. -□3 months, -20 to -70 °C under sterile conditions after reconstitution.

Images



Note: This product is for in vitro research use only and is not intended for use in humans or animals.