

# Recombinant Rat Interleukin-9

Catalog No: #AP60233

Package Size: #AP60233-1 10ug #AP60233-2 100ug #AP60233-3 500ug


Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Rat Interleukin-9
Host Species	Escherichia coli.
Purification	> 97 % by SDS-PAGE and HPLC analyses.
Calculated MW	Approximately 14.3 kDa, a single non-glycosylated polypeptide chain containing 127 amino acids.
Target Sequence	MQRCSTSWG I QHTSYLIENL KDDPSSK CSC SANVT SCLCL PIPSD DCTTP CFQEGMSQVT NATQQSKFSP FFFRVKRIVE TLKSNKCQFF SCEKPCNQTT AGNTVSFLKS LLKTFQKTEV QVQRSRA
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. <ul style="list-style-type: none"> <li>-□ A minimum of 12 months from date of receipt, when stored at &amp;le;-20 °C as supplied.</li> <li>-□ 1 month, 2 to 8 °C under sterile conditions after reconstitution.</li> <li>-□ 3 months, -20 to -70 °C under sterile conditions after reconstitution.</li> </ul>

## Images



Note: This product is for in vitro research use only and is not intended for use in humans or animals.