

## Recombinant Rat Platelet-derived Growth Factor-AA

Catalog No: #AP60256




Package Size: #AP60256-1 10ug #AP60256-2 100ug #AP60256-3 500ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Rat Platelet-derived Growth Factor-AA
Host Species	Escherichia coli.
Purification	> 98 % by SDS-PAGE and HPLC analyses.
Calculated MW	Approximately 25.3 kDa, a disulfide-linked homodimeric protein containing two 111 amino acid residues polypeptide (A chain).
Target Sequence	MSIEEAIPAV CKTRTVIYEI PRSQVDPTSA NFLIWPPCVE VKRCTGCCNT SSVKQCPSRV HHRSVKVAKV EYVRKKPKLK EVQVRLEEHL ECACATSNLN PDHREEETDV R
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in 30 % Acetonitrile and 0.1 % TFA.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. <input type="checkbox"/> A minimum of 12 months from date of receipt, when stored at &le;-20 °C as supplied. <input type="checkbox"/> 1 month, 2 to 8 °C under sterile conditions after reconstitution. <input type="checkbox"/> 3 months, -20 to -70 °C under sterile conditions after reconstitution.

## Images



Note: This product is for in vitro research use only and is not intended for use in humans or animals.