

Recombinant Rat Midkine

Catalog No: #AP60262

Package Size: #AP60262-1 5ug #AP60262-2 100ug #AP60262-3 500ug


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Rat Midkine
Host Species	Escherichia coli.
Purification	> 97 % by SDS-PAGE and HPLC analyses.
Calculated MW	Approximately 13.2 kDa, a single non-glycosylated polypeptide chain containing 120 amino acids.
Target Sequence	VAKKKDKVKK GSECSEWTWG PCTPSSKDCG MGFREGTCGA QTQRIHCKVP CNWKKEFGAD CKYKFESWGA CDGSTGTKAR QGTLKKARYN AQCQETIRVT KPCTSKTKSK AKAKKGKGD
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH7.4.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. -□ A minimum of 12 months from date of receipt, when stored at ≤-20 °C as supplied. -□ 1 month, 2 to 8 °C under sterile conditions after reconstitution. -□ 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Images



Note: This product is for in vitro research use only and is not intended for use in humans or animals.