

Recombinant Porcine Interleukin-1 Receptor Antagonist Protein

Catalog No: #AP60270


Package Size: #AP60270-1 5ug #AP60270-2 100ug #AP60270-3 500ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Porcine Interleukin-1 Receptor Antagonist Protein
Host Species	Escherichia coli.
Purification	> 97 % by SDS-PAGE and HPLC analyses.
Calculated MW	Approximately 17.1 kDa, a single non-glycosylated polypeptide chain containing 152 amino acids.
Target Sequence	HPLGKRPCR M QAFRIWDVNQ KTFYLRNNQL VAGYLQGPNT KLEEKIDVVP VEPHFVFLGI HGGKLCCLSCV KSGDEMKLQL DAVNITDLRK NSEQDKRFTF IRSDSGPTTS FESAACPGWF LCTALEADQP VGLTNTPKAA VKVTKFYFQQ DQ
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in 1 x PBS, pH 7.4, 1 mM DTT.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. <ul style="list-style-type: none"> -□ A minimum of 12 months from date of receipt, when stored at &le;-20 °C as supplied. -□ 1 month, 2 to 8 °C under sterile conditions after reconstitution. -□ 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Images



Note: This product is for in vitro research use only and is not intended for use in humans or animals.