

## Recombinant Equine Interleukin-2 Cys141Ser

Catalog No: #AP60275


Package Size: #AP60275-1 10ug #AP60275-2 100ug #AP60275-3 500ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Equine Interleukin-2 Cys141Ser
Host Species	Escherichia coli.
Purification	> 97 % by SDS-PAGE and HPLC analyses.
Calculated MW	Approximately 14.9 kDa, a single non-glycosylated polypeptide chain containing 129 amino acids.
Target Sequence	APTSSSKRET QQQLKQLQMD LKLLLEGVNN NKNPKLSKML TFKINMPKKA TELKHLQCLE EELKPLEEML KNFLSKDIKE LMSNINVTVL GLKGSETRFT CEYDDEGTI VEFLNKWITF SQSIFSTMT
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 6.5.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. -□ A minimum of 12 months from date of receipt, when stored at &le;-20 °C as supplied. -□ 1 month, 2 to 8 °C under sterile conditions after reconstitution. -□ 3 months, -20 to -70 °C under sterile conditions after reconstitution.

## Images



Note: This product is for in vitro research use only and is not intended for use in humans or animals.