

Recombinant Human Stromal-Cell Derived Factor-1 gamma/CXCL12 gamma

Catalog No: #AP60290


Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP60290-1 10ug #AP60290-2 100ug #AP60290-3 500ug

Description

Product Name	Recombinant Human Stromal-Cell Derived Factor-1 gamma/CXCL12 gamma
Host Species	Escherichia coli.
Purification	> 96 % by SDS-PAGE and HPLC analyses.
Calculated MW	Approximately 11.6 kDa, a single non-glycosylated polypeptide chain containing 99 amino acid residues.
Target Sequence	GKPVLSYRC PCRFFESHVA RANVKHLKIL NTPNCALQIV ARLKNNNRQV CIDPKLKWIQ EYLEKALNKG RREEKVGKKE KIGKKKRQKK RCAAQKRKN
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. <ul style="list-style-type: none"> -□ A minimum of 12 months from date of receipt, when stored at &le;-20 °C as supplied. -□ 1 month, 2 to 8 °C under sterile conditions after reconstitution. -□ 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Images



Note: This product is for in vitro research use only and is not intended for use in humans or animals.