

Recombinant Human VEGF Co-regulated Chemokine 1/CXCL17

Catalog No: #AP60293


Package Size: #AP60293-1 5ug #AP60293-2 100ug #AP60293-3 500ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human VEGF Co-regulated Chemokine 1/CXCL17
Host Species	Escherichia coli.
Purification	> 96 % by SDS-PAGE and HPLC analyses.
Calculated MW	Approximately 11.5 kDa, a single non-glycosylated polypeptide chain containing 98 amino acids.
Target Sequence	SSLNPGVARG HRDRGQASRR WLQEGGQECE CKDWFLRAPR RKFMTVSGLP KKQCPCDHFK GNVKKTRHQR HHRKPNKHSR ACQQFLKQCQ LRSFALPL
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4, with 3 % Trehalose.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. <input type="checkbox"/> A minimum of 12 months from date of receipt, when stored at ≤-20 °C as supplied. <input type="checkbox"/> 1 month, 2 to 8 °C under sterile conditions after reconstitution. <input type="checkbox"/> 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Images



Note: This product is for in vitro research use only and is not intended for use in humans or animals.