

# Recombinant Murine Macrophage Inflammatory Protein-1 alpha/CCL3

Catalog No: #AP60345


Package Size: #AP60345-1 10ug #AP60345-2 100ug #AP60345-3 500ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Murine Macrophage Inflammatory Protein-1 alpha/CCL3
Host Species	Escherichia coli.
Purification	> 98 % by SDS-PAGE and HPLC analyses.
Calculated MW	Approximately 7.9 kDa protein containing 69 amino acid residues, including the four highly conserved cysteine residues present in CC chemokines.
Target Sequence	APYGADTPTA CCFSYSRKIP RQFIVDYFET SSLCSQPGVI FLTKRNRQIC ADSKETWVQE YITDLELNA
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in 30 % Acetonitrile and 0.1 % TFA.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. <ul style="list-style-type: none"> <li>-□ A minimum of 12 months from date of receipt, when stored at &amp;le;-20 °C as supplied.</li> <li>-□ 1 month, 2 to 8 °C under sterile conditions after reconstitution.</li> <li>-□ 3 months, -20 to -70 °C under sterile conditions after reconstitution.</li> </ul>

## Images



Note: This product is for in vitro research use only and is not intended for use in humans or animals.