

Recombinant Human Parathyroid Hormone 7-84

Catalog No: #AP60398




Package Size: #AP60398-1 10ug #AP60398-2 100ug #AP60398-3 500ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Parathyroid Hormone 7-84
Host Species	E.coli
Purification	> 98 % by SDS-PAGE and HPLC analyses.
Calculated MW	Approximately 8.8 kDa, a single non-glycosylated polypeptide chain containing 78 amino acids.
Target Sequence	LMHNLGKHLN SMERVEWLRK KLQDVHNFVA LGAPLAPRDA GSQRPRKKED NVLVESHEKS LGEADKADVN VLTKAKSQ
Formulation	LyophilizedB fromB aB 0.2B umB filteredB concentratedB solutionB inB PBS.B
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. -□ A minimum of 12 months from date of receipt, when stored at ≤-20 °C as supplied. -□ 1 month, 2 to 8 °C under sterile conditions after reconstitution. -□ 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Images



Note: This product is for in vitro research use only and is not intended for use in humans or animals.