

# Recombinant Human Ubiquitin-conjugating Enzyme E2 M

Catalog No: #AP60439

Package Size: #AP60439-1 5ug #AP60439-2 100ug #AP60439-3 500ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Ubiquitin-conjugating Enzyme E2 M
Host Species	Escherichia coli.
Purification	> 95 % by SDS-PAGE and HPLC analyses.
Calculated MW	Approximately 20.9 kDa, a single non-glycosylated polypeptide chain containing 183 amino acids.
Target Sequence	MIKLFSLKQQ KKEEESAGGT KGSSKKASAA QLRIQKDINE LNLPKTCDIS FSDPDDLLNF KLVICPDEGF YKSGKFVFSF KVGQGYPHDP PKVKCETMVY HPNIDLEGNV CLNILREDWK PVLTIINSIY GLQYLFLEPN PEDPLNKEAA EVLQNNRRLF EQNVQRSMRG GYIGSTYFER CLK
Formulation	A 0.2 µm filtered concentrated solution in 50 mM HEPES, pH 6.5, with 125 mM NaCl, 10 % Glycerol, 5 % Trehalose, 1 mM DTT.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. <ul style="list-style-type: none"> <li>-□ A minimum of 12 months from date of receipt, when stored at &amp;le;-20 °C as supplied.</li> <li>-□ 1 month, 2 to 8 °C under sterile conditions after reconstitution.</li> <li>-□ 3 months, -20 to -70 °C under sterile conditions after reconstitution.</li> </ul>

Note: This product is for in vitro research use only and is not intended for use in humans or animals.