

# Recombinant Human Cyclin-Dependent Kinase Inhibitor 2A, Isoform 1-TAT

Catalog No: #AP60462


Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP60462-1 5ug #AP60462-2 100ug #AP60462-3 500ug

## Description

Product Name	Recombinant Human Cyclin-Dependent Kinase Inhibitor 2A, Isoform 1-TAT
Host Species	Escherichia coli.
Purification	> 95 % by SDS-PAGE and HPLC analyses.
Other Names	CDK4I, MTS-1, p16INK4A
Calculated MW	Approximately 18.0 kDa, a single non-glycosylated polypeptide chain containing 167 amino acids.
Target Sequence	EPAAGSSMEP SADWLATAAA RGRVEEVRAL LEAGALPNAP NSYGRRPIQV MMMGSARVAE LLLLHGAEPN CADPATLTRP VHDAAREGFL DTLVVLHRAG ARLDVRDAWG RLPVDLAEEL GHRDVARYLR AAAGGTRGSN HARIDAAEGP SDIPDGYGRK KRRQRRR
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in 2 x PBS, pH 7.0.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. -□ A minimum of 12 months from date of receipt, when stored at &le;-20 °C as supplied. -□ 1 month, 2 to 8 °C under sterile conditions after reconstitution. -□ 3 months, -20 to -70 °C under sterile conditions after reconstitution.

## Images



## Background

Cyclin-dependent kinase inhibitors (CDKIs) are proteins that bind to and inhibit the activity of CDKs. Two major classes of CDK inhibitors have been identified. The p16 family (p15, p16, p18 and p19) binds to and inhibits the activities of CDK4 and CDK6. The p21 family (p21, p27, p28 and p57) can bind to broad range of CDK-cyclin complexes and inhibit their activities. CDKIs are capable of suppressing growth, and several lines of evidence strongly suggest that at least some CDKIs may be tumor suppressor proteins.

p16-INK4A is the member of p16 family and is encoded by CDKN2A gene in humans. It has three isoforms, which are widely expressed but not detected in brain or skeletal muscle, except that isoform 3 is pancreas-specific. Defects in p16INK4A are a cause of Li-Fraumeni syndrome (LFS) and melanoma-astrocytoma syndrome (MASTS).

The TAT transduction peptide can help rHuP16-INK4a with crossing all kind of biomembranes in vivo and has no effect on the super structure of the protein.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.