

Recombinant Human Epididymis Protein 4

Catalog No: #AP60484



Package Size: #AP60484-1 10ug #AP60484-2 100ug #AP60484-3 500ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Epididymis Protein 4
Host Species	Escherichia coli.
Purification	> 95 % by SDS-PAGE.
Calculated MW	Approximately 10.0 kDa, a single polypeptide chain containing 94 amino acids. But it migrates with an apparent molecular mass of 16.9 kDa in SDS-PAGE.
Target Sequence	EKTGVCPQLQ ADQNCTQECV SDSECADNLK CCSAGCATFC SLPNDKEGSC PQVNINFPQL GLCRDQCQVD SQCPGQMKCC RINGCGKVSCV TPNF
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.0.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. -□ A minimum of 12 months from date of receipt, when stored at ≤-20 °C as supplied. -□ 1 month, 2 to 8 °C under sterile conditions after reconstitution. -□ 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.