

# Recombinant Murine Desert Hedgehog Homolog Cys231Ilelle

Catalog No: #AP60491


Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP60491-1 5ug #AP60491-2 100ug #AP60491-3 500ug

## Description

Product Name	Recombinant Murine Desert Hedgehog Homolog Cys231Ilelle
Host Species	Escherichia coli.
Purification	> 98 % by SDS-PAGE and HPLC analyses.
Calculated MW	Approximately 20.0 kDa, a single non-glycosylated polypeptide chain containing 177 amino acids.
Target Sequence	IIGPGRGPVG RRRYVRKQLV PLYKQFVPS MPERTLGASG PAEGRVTRGS ERFRDLVPNY NPDIIIFKDEE NSGADRLMTE RCKERVNALA IAVMNMWPGV RLRVTEGWDE DGHHAQDSLH YEGRALDITT SDRDRNKYGL LARLAVEAGF DWVYYESRNH IHVSVKADNS LAVRAGG
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4, 1 mM DTT, 0.05 % Tween-80.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. <ul style="list-style-type: none"> <li>-□ A minimum of 12 months from date of receipt, when stored at &amp;le;-20 °C as supplied.</li> <li>-□ 1 month, 2 to 8 °C under sterile conditions after reconstitution.</li> <li>-□ 3 months, -20 to -70 °C under sterile conditions after reconstitution.</li> </ul>

## Images



Note: This product is for in vitro research use only and is not intended for use in humans or animals.