

# Recombinant Hirudin

Catalog No: #AP60500



Package Size: #AP60500-1 10ug #AP60500-2 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Hirudin
Host Species	Escherichia coli
Purification	> 96 % HPLC analyses.
Target Sequence	LTVVYTDCTESGQNLCLCEGSNVCGQGKNCILGSDGEKNQCVTGEGTPGPQSHNDGDFEEP EEYL
Formulation	Lyophilized from a 0.2 µm filtered solution of 20 mM PB, pH 7.0, containing 2 % mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. -□ A minimum of 12 months from date of receipt, when stored at &le;-20 °C as supplied. -□ 1 month, 2 to 8 °C under sterile conditions after reconstitution. -□ 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.