

Recombinant Single-stranded DNA Binding Protein

Catalog No: #AP60509




Package Size: #AP60509-1 100ug #AP60509-2 500ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Single-stranded DNA Binding Protein
Host Species	Escherichia coli
Purification	> 95 % by SDS-PAGE and HPLC analyses.
Calculated MW	Approximately 16.1 kDa, a single non-glycosylated polypeptide chain containing 148 amino acids.
Target Sequence	MEEKVGNLKP NMESVNVTVR VLEASERQI QTKNGVRTIS EAIVGDETGR VKLTLWGKHA GSIKEGQVVK IENAWTTAFK GQVQLNAGSK TKIAEASEDG FPESQIPEN TPTAPQQMRG GGRGFRGGGR RYGRRGGRRQ ENEEGEEE
Formulation	0.2 µm filtered concentrated solution in 20 mM Tris, pH 7.4, 200 mM NaCl, 50 % Glycerol, 1 mM EDTA, 0.5 mM DTT.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. <input type="checkbox"/> A minimum of 12 months from date of receipt, when stored at ≤-20 °C as supplied. <input type="checkbox"/> 1 month, 2 to 8 °C under sterile conditions after reconstitution. <input type="checkbox"/> 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Images



Note: This product is for in vitro research use only and is not intended for use in humans or animals.