



Recombinant Mouse Long-chain-fatty-acid--CoA
ligase 6(Acsl6),partial

Catalog No: #AP89657

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Package Size: #AP89657-1 20ug #AP89657-2 100ug #AP89657-3 1mg

Description

Product Name	Recombinant Mouse Long-chain-fatty-acid--CoA ligase 6(Acsl6),partial
Host Species	E.coli
Purification	>85% (SDS-PAGE)
Target Name	His-tagged
Other Names	Acsl6; Fac16Long-chain-fatty-acid--CoA ligase 6; EC 6.2.1.3; Arachidonate--CoA ligase; EC 6.2.1.15; Long-chain acyl-CoA synthetase 6; LACS 6
Target Species	Mus musculus (Mouse)
Formulation	Protein Length:Partial Expression Region: 46-697aa HRPKALQPPCNLLKQSEEVEDGGGARRSVIGGCTQLLTHYYDDARTMYQVFRRGLSISGNGPCLGFRKPEQP YQWLSYQEVAKRAEFLGSGLLQHDKVGTQFVGVFAQNRPEWIIAELACYTYSMVVPLYDTLPGSISYIIN TADICTVIVDKPHKATLLLEHVERKETPGLKLVILMEPFEDALRERGGKCGVDIKSMQAIEDCGRENHHAPVPP RPDDLIVCFTSGTTGNPKGAMLTHGNVADFSGFLKVTEKVIFPRQDDVLISFLPLAHMFERVIQSVVYCHGG RVGFFQGDIRLLSDDMKALRPTIFPVPRLLNRMYSKIFHQADTSLKRWLLEFAAKRKQAEVRSGIIRNNSIWD ELFFNKIQASLGGHVRMIVTGAAPASPTVLGFLRAALGCQVYEGYGQTECTAGCTFTTPGDWTSGHVGAPLP CNHIKLV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.