

Fish Tumor Necrosis Factor Alpha(TNF α)ELISA Kit

Catalog No: #EK20560



Package Size: #EK20560-1 48T #EK20560-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Fish Tumor Necrosis Factor Alpha(TNF α)ELISA Kit
Applications	ELISA
Species Reactivity	Fish

Application Details

Detect Range:7.5pg/ml-90pg/ml

Sensitivity:1.2pg/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.