

Myc tag mouse Monoclonal Antibody APC Conjugated

Catalog No: #C92847APC



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Myc tag mouse Monoclonal Antibody APC Conjugated
Host Species	Mouse
Clonality	Monoclonal
Purification	affinity purified by Protein G
Applications	ICC IF
Species Reactivity	Hu Myc tag
Conjugates	APC
Target Name	Myc tag
Other Names	Myc tag
Accession No.	Uniprot:P01106Gene ID:4609
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.