

# Human THO complex subunit 4 (THOC4) ELISA Kit

Catalog No: #EK20582



Package Size: #EK20582-1 48T #EK20582-2 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human THO complex subunit 4 (THOC4) ELISA Kit
Applications	ELISA
Species Reactivity	Human

## Application Details

Detect Range:0.312-20 ng/mL
Sensitivity:0.125ng/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 &mu;L
Assay Time:1-4.5h
Detection wavelength:450 nm
Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.