

## TRIF Rabbit mAb

Catalog No: #57354

Package Size: #57354-1 50ul #57354-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	TRIF Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human TRIF
Other Names	TRIF; IIAE6; MyD88-3; PRVTIRB; TICAM-1
Accession No.	Swiss-Prot:Q8IUC6 Gene ID:148022
Calculated MW	76 kDa
SDS-PAGE MW	98 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


## Application Details

IP:1/20WB: 1/500-1/1000

## Images



Western blot analysis of TRIF in rat testis lysates using TRIF antibody.



Western blot analysis of TRIF in K562, 3T3, HeLa lysates using TRIF antibody.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.