

MRPS18B Rabbit mAb

Catalog No: #57382

Package Size: #57382-1 50ul #57382-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	MRPS18B Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP ICC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human MRPS18B
Other Names	PTD017; S18amt; C6orf14; HSPC183; MRPS18-2; HumanS18a; MRP-S18-2
Accession No.	Swiss-Prot:Q9Y676 Gene ID:28973
Calculated MW	29 kDa
SDS-PAGE MW	29 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IP:1/20WB: 1/500-1/1000ICC: 1/50-1/200

Images



Western blot analysis of MRPS18B in Raw264.7 lysates using MRPS18B antibody.



Western blot analysis of MRPS18B in A549, HL-60, U251 lysates using MRPS18B antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.