

Mouse Ifngr1/ Interferon gamma receptor 1 ELISA Kit

Catalog No: #EK21127

Package Size: #EK21127 96T



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Ifngr1/ Interferon gamma receptor 1 ELISA Kit
Applications	ELISA
Species Reactivity	Mouse
Other Names	Interferon gamma receptor alpha-chain, IFN-gamma-R-alpha, CD119, IFN-gamma receptor 1, Ifngr

Application Details

Detect Range:0.312-20ng/mL
Sensitivity:0.156ng/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm
Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.