

Rat Alpha-2-Plasmin Inhibitor (a2PI) ELISA Kit

Catalog No: #EK20575



Package Size: #EK20575-1 48T #EK20575-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Alpha-2-Plasmin Inhibitor (a2PI) ELISA Kit
Applications	ELISA
Species Reactivity	Rat
Other Names	SERPINF2

Application Details

Detect Range:0.312-20ng/mL

Sensitivity:0.115 ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.