

Human Paired immunoglobulin-like type 2 receptor alpha (PILRA) ELISA Kit

Catalog No: #EK20576



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #EK20576-1 48T #EK20576-2 96T

Description

Product Name	Human Paired immunoglobulin-like type 2 receptor alpha (PILRA) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:78.1-5000pg/mL

Sensitivity:27pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.