

Recombinant Human Charged multivesicular body protein 1a(CHMP1A),Escherichia Coli

Catalog No: #AP89666

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP89666-1 20ug #AP89666-2 100ug #AP89666-3 1mg

Description

Product Name	Recombinant Human Charged multivesicular body protein 1a(CHMP1A),Escherichia Coli
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	>85% (SDS-PAGE)
Immunogen Description	Expression Region:1-196Sequence Info:full length protein
Target Name	His-tagged
Other Names	Charged multivesicular body protein 1a; CHM1A_HUMAN; CHMP1a; Chromatin modifying protein 1a; Chromatin-modifying protein 1a; hVps46 1; hVps46-1; KIAA0047; Metalloprotease 1; PCOLN3; Procollagen (type III) N endopeptidase; PRSM1; Vacuolar protein sorting 46 1; Vacuolar protein sorting-associated protein 46-1; Vps46 1; Vps46-1
Target Species	Homo sapiens (Human)
Tag Info	N-terminal His-tagged Tag-Free
Target Sequence	MDDTLFQLKF TAKQLEKLAK KAEKDSKAEQ AKVKKALLQK NVECARVYAE NAIRKKNQGV NWLRMASRVD AVASKVQTAV TMKGVTKNMA QVTKALDKAL STMDLQKVSS VMDFEQVQV NLDVHTSVME DSMSSATTLT TPQEQVDSL I MQIAEENGLE VLDQLSQLPE GASAVGESSV RSQEDQLSRR LAALRN
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.