

Recombinant Human Integrin alpha-M (ITGAM), partial

Catalog No: #AP89675

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP89675-1 20ug #AP89675-2 100ug #AP89675-3 1mg

Description

Product Name	Recombinant Human Integrin alpha-M (ITGAM), partial
Brief Description	Recombinant Human Integrin alpha-M (ITGAM), partial
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Target Name	His-tagged
Target Species	Homo sapiens (Human)
Calculated MW	57.1 kDa
Tag Info	N-terminal 6xHis-GST-tagged
Target Sequence	"Protein Length:Partial Expression Region: 46-150aa VLFQGPLGSPEFRTVVVGAPQEIVAANQRGSLYQCDYSTGSCEPIRLQVPVEAVNMSLGLSLAATTSPQLLA CGPTVHQTCSENTYVKGLCFLFGSNLRQQPQKFPEALRGCPQEDSDGRTRAPPPPLRSGCQSPKGSVGC CHRAITSITPWGLTGLEGFFAERRNYIRIGEWDAPCSGALSAAGVVVTRSVTATLASALAPAPFAFFPSFLATFA GFPRQALNRGLPLGFRFSALRHL"
Storage	"The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C."

Note: This product is for in vitro research use only and is not intended for use in humans or animals.