

Rhodopsin Polyclonal Antibody

Catalog No: #48914



Package Size: #48914-1 50ul #48914-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rhodopsin Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;IHC;IF;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	Rhodopsin Polyclonal Antibody detects endogenous levels of Rhodopsin protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human Rhodopsin.
Calculated MW	42kD
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C 1 year

Application Details


WB 1:500 - 1:2000.

IHC 1:100 - 1:300.


ELISA: 1:5000..

IF 1:50-200

Images



Western Blot analysis of various cells using Rhodopsin Polyclonal Antibody



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using Rhodopsin Antibody. The picture on the right is blocked with the synthesized peptide.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.