

TRAP Mouse mAb

Catalog No: #60450

Package Size: #60450-1 50ul #60450-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	TRAP Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1, Kappa
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Applications	IHC;IF
Species Reactivity	Human;Mouse;Rat
Specificity	The antibody can specifically recognize human TRAP protein.
Immunogen Description	Synthesized peptide derived from human TRAP
Target Name	ACP5
Other Names	Tartrate-resistant acid phosphatase type 5 (TR-AP;EC 3.1.3.2;Tartrate-resistant acid ATPase;TrATPase;Type 5 acid phosphatase)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


IHC-p 1:200-400,

IF 1:200-400

Images



Human lung tissue was stained with Anti-TRAP Antibody



Human giant cell tumor tissue was stained with Anti-TRAP Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.