

# Mouse Pf4/ Platelet factor 4 ELISA Kit

Catalog No: #EK21120

Package Size: #EK21120-1 96T



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Pf4/ Platelet factor 4 ELISA Kit
Brief Description	Mouse Pf4/ Platelet factor 4 ELISA Kit
Applications	ELISA
Species Reactivity	Mouse
Other Names	C-X-C motif chemokine 4, PF-4, Cxcl4, Scyb4

## Application Details

Detect Range:78-5000pg/mL

Sensitivity:39.8pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.