

Porcine B-Cell Leukemia/Lymphoma 2 (BCL2) ELISA Kit

Catalog No: #EK20610

Package Size: #EK20610-1 48T #EK20610-2 96T



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Porcine B-Cell Leukemia/Lymphoma 2 (BCL2) ELISA Kit
Brief Description	Porcine B-Cell Leukemia/Lymphoma 2 (BCL2) ELISA Kit
Applications	ELISA
Species Reactivity	Porcine

Application Details

Detect Range:31.2-2000 pg/mL

Sensitivity:13.7pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.