

Porcine Caspase 9 (CASP9) ELISA Kit

Catalog No: #EK20614



Package Size: #EK20614-1 48T #EK20614-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|-------------------------------------|
| Product Name | Porcine Caspase 9 (CASP9) ELISA Kit |
| Brief Description | Porcine Caspase 9 (CASP9) ELISA Kit |
| Applications | ELISA |
| Species Reactivity | Porcine |

Application Details

Detect Range:0.156-10 ng/mL

Sensitivity:0.057ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.