

Mouse Collagen Type IX Alpha 1 (COL9a1) ELISA Kit

Catalog No: #EK20618



Package Size: #EK20618-1 48T #EK20618-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Collagen Type IX Alpha 1 (COL9a1) ELISA Kit
Brief Description	Mouse Collagen Type IX Alpha 1 (COL9a1) ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:15.6-1000pg/mL

Sensitivity:6.7pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.