

CXCR2 Conjugated Antibody

Catalog No: #C29591



Package Size: #C29591-AF350 100ul #C29591-AF405 100ul #C29591-AF488 100ul #C29591-AF555 100ul #C29591-AF594 100ul #C29591-AF647 100ul #C29591-AF680 100ul #C29591-AF750 100ul #C29591-Biotin 100ul #C29591-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	CXCR2 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human CXCR2.
Other Names	CXCR2;CD182;CDw128b;CMKAR2;IL8R2;IL8RA;IL8RB
Calculated MW	40kDa
SDS-PAGE MW	41kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at 4Λ C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.