

Rat Fibrinogen (Fg) ELISA Kit

Catalog No: #EK20620

Package Size: #EK20620-1 48T #EK20620-2 96T



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|-------------------------------|
| Product Name | Rat Fibrinogen (Fg) ELISA Kit |
| Brief Description | ELISA Kit |
| Applications | ELISA |
| Species Reactivity | Rat |

Application Details

Detect Range:78.1-5000 ng/mL

Sensitivity:29ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.