

Rat Nuclear factor of activated T-cells, cytoplasmic 1 (NFATC1) ELISA Kit

Catalog No: #EK20635

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #EK20635-1 48T #EK20635-2 96T

Description

Product Name	Rat Nuclear factor of activated T-cells, cytoplasmic 1 (NFATC1) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Rat

Application Details

Detect Range:0.156-10ng/mL
Sensitivity:0.064ng/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.