

Recombinant Mouse L-amino-acid oxidase (Il4i1)

Catalog No: #AP89674



Package Size: #AP89674-1 20ug #AP89674-2 100ug #AP89674-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse L-amino-acid oxidase (Il4i1)
Brief Description	Recombinant Proteins
Host Species	E.coli
Purification	>85% (SDS-PAGE)
Immunogen Description	Expression Region: 22-630Sequence Info:Full Length of Mature Protein
Target Name	His-tagged
Other Names	Il4i1; Fig1; L-amino-acid oxidase; LAAO; LAO; EC 1.4.3.2; Interleukin-4-induced protein 1; IL4-induced protein 1; Protein Fig-1; mFIG1
Target Species	Mus musculus (Mouse)
Tag Info	N-terminal His-tagged Tag-Free
Target Sequence	DWKAASSLN PIEKCMEDHD YEQLLKVVTL GLNRTSKPQK VVVVGAGVAG LVAAKMLSDA GHKVTILEAD NRIGGRITF RDEKTGWIGE LGAMRMPSSH RILHKLCRTL GLNLTQFTQY DENTWTEVHN VKLRNYVVEK MPEKLGYNLN NRERGHSPED IYQMALNKAF KDLKALGCKK AMNKFNKHTL LEYLLEEGNL SRPAVQLLGD VMSEEGFFYL SFAEALRAHA CLSDRLRYSR IVGGWDL LPR ALLSSLSGAL LLNAPVVSIT QGRNDVRVHI ATSLHSEKTL TADVLLTAS GPALQRITFS PPLTRKRQEA LRALHYVAAS KVFLSFRFPF WHEEHIEGGH SNTDRPSRLI FYPARGEGSL LLASYTWSDA AAPFAGLSTD QTLRLVLQDV AALHGPVVFR LWDGRGVVKR WAEDPHSQGG FVVQPPLYGR EAEDYDWSAP FGRIYFAGEH TALPHGWVET AVKSGLRAAV RINNNGYGE VDPQMMEHAY AEANYLDQYP EGERPEEQQA REEVSPDEQE PSHKLLVET SPEGQQHAFV EAIPELQGHV FVETVPQKEG HAHQNIYPSE HVQVHGEVIP EWHGHGGSST PQMHRVGDHS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.