

## CD204 Polyclonal Antibody

Catalog No: #48940

Package Size: #48940-1 50ul #48940-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CD204 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC;IF;ELISA
Species Reactivity	Human;Mouse;Rat
Immunogen Description	The antiserum was produced against synthesized peptide derived from the Internal region of human MSR1.
Target Name	MSR1 SCARA1
Other Names	Macrophage scavenger receptor types I and II (Macrophage acetylated LDL receptor I and II;Scavenger receptor class A member 1;CD antigen CD204)
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

IHC-p 1:50-200,

ELISA 1:10000-20000.

IF 1:50-200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.