

Rat Eukaryotic translation initiation factor 4E-binding protein 1 (EIF4EBP1) ELISA Kit

Catalog No: #EK20633



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #EK20633-1 48T #EK20633-2 96T

Description

Product Name	Rat Eukaryotic translation initiation factor 4E-binding protein 1 (EIF4EBP1) ELISA Kit
Brief Description	Rat Eukaryotic translation initiation factor 4E-binding protein 1 (EIF4EBP1) ELISA Kit
Applications	ELISA
Species Reactivity	Rat

Application Details

Detect Range:0.156-10ng/mL

Sensitivity:0.062ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.