

IL-2R gamma Polyclonal Antibody

Catalog No: #C91756



Package Size: #C91756-1 50ul #C91756-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	IL-2R gamma Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt, ICC
Species Reactivity	Hu Ms Rat
Conjugates	Unconjugated
Target Name	IL-2R gamma
Other Names	IL-2R gamma
Accession No.	Uniprot:P31785Gene ID:3561
Calculated MW	41
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt $1 \mu\text{g} / \text{test}$

ICC $1:100$

Note: This product is for in vitro research use only and is not intended for use in humans or animals.