

HumanB 8-isoB prostaglandinB F2aB (8-iso-PGF2a)B ELISAB Kit

Catalog No: #EK20643



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #EK20643-1 48T #EK20643-2 96T

Description

Product Name	HumanB 8-isoB prostaglandinB F2aB (8-iso-PGF2a)B ELISAB Kit
Brief Description	HumanB 8-isoB prostaglandinB F2aB (8-iso-PGF2a)B ELISAB Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:125pg/mL-8000 pg/mL

Sensitivity:46pg/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Detection Method:competitive inhibition

Note: This product is for in vitro research use only and is not intended for use in humans or animals.