

Human 3-oxoacyl-[acyl-carrier-protein] synthase, mitochondrial (OXSM) ELISA Kit

Catalog No: #EK20637



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #EK20637-1 48T #EK20637-2 96T

Description

Product Name	Human 3-oxoacyl-[acyl-carrier-protein] synthase, mitochondrial (OXSM) ELISA Kit
Brief Description	Human 3-oxoacyl-[acyl-carrier-protein] synthase, mitochondrial (OXSM) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.156-10ng/mL

Sensitivity:0.067/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.