

# Mouse Neural Cell Adhesion Molecule (NCAM) ELISA Kit

Catalog No: #EK20660

Package Size: #EK20642-1 #EK20642-2



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Mouse Neural Cell Adhesion Molecule (NCAM) ELISA Kit |
| Brief Description  | Mouse Neural Cell Adhesion Molecule (NCAM) ELISA Kit |
| Applications       | ELISA  |
| Species Reactivity | Mouse  |

## Application Details

|  |
|--|
| Detect Range: 78.1-5000pg/mL   |
| Sensitivity:29pg/mL  |
| Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| Sample Volume: 50 &mu;L  |
| Assay Time:1-4.5h  |
| Detection wavelength:450 nm  |
| Detection Method:sandwich  |

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.