

Recombinant Human NAD-dependent protein deacetylase sirtuin-3, mitochondrial (SIRT3), partial

Catalog No: #AP89676

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP89676-1 20ug #AP89676-2 100ug #AP89676-3 1mg

Description

Product Name	Recombinant Human NAD-dependent protein deacetylase sirtuin-3, mitochondrial (SIRT3), partial
Brief Description	Recombinant Human NAD-dependent protein deacetylase sirtuin-3, mitochondrial (SIRT3), partial
Host Species	E.coli
Purification	Greater than 95% as determined by SDS-PAGE.
Target Name	His-tagged
Other Names	(hSIRT3)(Regulatory protein SIR2 homolog 3)(SIR2-like protein 3)
Target Species	Homo sapiens (Human)
Calculated MW	32.3 kDa
Tag Info	C-terminal 6xHis-tagged
Target Sequence	"Protein Length:Partial Expression Region: 118-399aa SDKGKLSLQDVAELIRARACQRVVMVGAGISTPSGIPDFRSPGSGLYSNLQQYDLPYPEAIFELPFFFHNPKP FFTLAKELYPGNYKPNVTHYFLRLLHDKGLLLRLYTQIDGLERVSGIPASKLVEAHGTFASATCTVCQRPFPG EDIRADVMADRVPRCPVCTGVVKPDIVFFGEPLPQRFLHVVDFPMADLLLILGTSLEVEPFASLTEAVRSSVP RLLINRDLVGPLAWHPRSRDVAQLGDVVHGVESLVELLGWTEEMRDLVQRETGKLDGPKD"
Storage	"The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C."

Note: This product is for in vitro research use only and is not intended for use in humans or animals.