

CHI3L1 Antibody

Catalog No: #C00652

Package Size: #C00652-1 50ul #C00652-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CHI3L1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB, IHC-P, IHC-F, Flow-Cyt, IF, ELISA
Species Reactivity	Human, Mouse, Rat
Immunogen Description	KLH conjugated synthetic peptide derived from human CHI3L1: 301-383/383
Conjugates	Unconjugated
Target Name	CHI3L1
Other Names	GP39; ASRT7; GP-39; YKL4; CGP-39; YKL-4; YYL-4; HC-gp39; HCGP-3P; hCGP-39; Chitinase-3-like protein 1; 39 kDa synovial protein; Cartilage glycoprotein 39; CHI3L1
Accession No.	Swiss-Prot#P36222NCBI Gene ID1116
Calculated MW	42kD
Concentration	1mg/ml
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB □ 1:500-2000

IHC-P □ 1:100-500

IHC-F □ 1:100-500

Flow-Cyt □ 1 μg/Test

IF □ 1:100-500

ELISA □ 1:5000-10000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.